

# Strategy for Training

To implement your hybrid model on the CUBIC cluster, we are combining three specific frameworks: the **DINOv3 ViT-B** (Vision), the **CheXzero/OpenAI CLIP** (Text), and the **PLIP** (Training Strategy).

Here are the definitive parameters for our implementation:

## 1. Architectural Parameters

Because the vision and text encoders come from different families, the "bridge" between them is the most critical parameter.

Parameter	Value	Description
Vision Backbone	DINOv3 ViT-B/16	Pretrained ViT-Base with 16x16 patch size.
Vision Output Dim	768	The raw embedding size from DINOv3 ViT-B.
Text Backbone	CLIP-ViT-B/32	The text encoder used in the CheXzero/OpenAI repo.
Text Output Dim	512	The latent space dimension used by CLIP/CheXzero.
Projection Head	Linear Layer	A layer mapping the Vision 768 ->512 to match the Text side.
Input Resolution	224 x 224	Standard input for both encoders; uses Bicubic interpolation.

## 2. Training Strategy Parameters (PLIP-Style)

We are shifting from the CheXzero "Epoch" approach to the PLIP "Step" approach.

Parameter	Value	Description
Total Training Steps	25,000	The total iteration count derived from the PLIP paper.
Batch Size	64	Strictly following the CheXzero batch size for memory stability.
Eval Interval	500 steps	Frequency of zero-shot validation on the external dataset.
Loss Function	InfoNCE	Symmetrical contrastive loss with learnable temperature.
Temperature	0.07	Initial logit scaling (often learnable in CLIP-based models).

## 3. Optimizer Hyperparameters (CheXzero-Style)

We are retaining the original CheXzero optimization settings to ensure the learning dynamics remain stable for medical images.

Parameter	Value	Description
Optimizer	SGD	Stochastic Gradient Descent.
Learning Rate	$10^{-4}$	The baseline rate from the CheXzero paper.

<b>Momentum</b>	<b>0.9</b>	Standard momentum for medical image convergence.
<b>Weight Decay</b>	<b>0</b>	CheXzero does not typically use weight decay in the final recipe.
<b>Text Context Length</b>	<b>512 tokens</b>	Uses the distilled context length for long medical reports.

#### 4. Data Preprocessing & Augmentation

<b>Parameter</b>	<b>Value</b>	<b>Description</b>
<b>Resize Strategy</b>	<b>320 -&gt; 224</b>	Images stored at 320px, then cropped/resized to 224px.
<b>Normalization Mean</b>	[0.481, 0.457, 0.408]	Specific to OpenAI CLIP weights.
<b>Normalization Std</b>	[0.268, 0.261, 0.275]	Specific to OpenAI CLIP weights.
<b>Augmentations</b>	Random Crop, Flip	Standard PLIP/CheXzero visual augmentations.